

In Vivo Adjustment Feature in the Masimo Radical-7™

OVERVIEW

The *In Vivo Adjustment* feature in the Radical-7 lets clinicians manually adjust the displayed value of one or more clinical parameters to match that of a corresponding laboratory reference for continuous trending. To remind clinicians that the feature is active, an “offset” value is displayed with the adjusted value.


The *In Vivo Adjustment* feature is available for the following clinical parameters:

- > Oxygen saturation (SpO₂)
- > Total hemoglobin (SpHb®)
- > Carboxyhemoglobin (SpCO®)
- > Methemoglobin (SpMet®)

The *In Vivo Adjustment* feature helps reduce individual patient bias that is expected when comparing a noninvasive measurement to a laboratory reference.


IN VIVO ADJUSTMENT OVERVIEW



For optimal performance, use the *In Vivo Adjustment* feature at the start of patient monitoring and during periods of uninterrupted patient monitoring. When the *In Vivo Adjustment* is set to Yes, the feature is active (turned on) and a positive or negative offset value appears, as shown in the following illustration. When the sensor or the patient cable is disconnected from the monitor, the offset value(s) will be reset to 0 (zero). The *In Vivo Adjustment* feature remains active until reset to No by the user.




(A) Displayed value of one or more adjusted clinical parameters. (B) Offset value appears and indicates that *In Vivo Adjustment* is active. A positive value means that the displayed value is increased (relative to the laboratory reference value) and a negative value means the displayed value is decreased (relative to the laboratory reference value).


STEPS TO USE THE IN VIVO ADJUSTMENT FEATURE



1. During patient monitoring, record the time and the value of the parameter to be adjusted.
2. At the same time, to verify the parameter, obtain a blood sample.
3. When the results return, subtract the recorded parameter value from the laboratory reference value. If the two values are identical, an adjustment is not needed.
4. On the monitor, go to the *In Vivo* screen:  > **In Vivo**


5. On the *In Vivo* screen, on the *In Vivo Adjustment* option, press  and use  or  to highlight **Yes**.





6. To confirm changes, select .
7. On the *In Vivo* screen, select the parameter you want to adjust, and select .



8. With the parameter highlighted, use  or  to enter the offset value. The offset value is the difference between the laboratory reference and the displayed parameter.



9. To confirm changes, select .
10. Repeat these steps as necessary for any of the parameters that are active.
11. When all adjustments are complete, press .

WARNING

CONFIRM OFFSET VALUE(S) PERIODICALLY AS THE DIFFERENCE BETWEEN THE DISPLAYED PARAMETER VALUE AND THE LABORATORY REFERENCE VALUE MAY VARY OVER TIME.

CAUTION

Do not use *In Vivo Adjustment* if the monitor displays a *Low SpHb SIQ* message.

MAXIMUM OFFSET FOR IN VIVO ADJUSTMENT

- > %SpO₂ (± 6% in increments of ± 1%)
- > SpCO (± 9% in increments of ± 1%)
- > SpMet (± 3% in increments of ± 0.1%)
- > SpHb (± 3 g/dL in increments of ± 0.1 g/dL)